

# Dental milling and grinding machines from vhf.

Find the perfect machine for  
your lab and practice.



■ EASE CLASS

■ PERFORMANCE CLASS

■ HIGH END CLASS

**vhf**  
CREATING PERFECTION

EASE CLASS

PERFORMANCE CLASS

HIGH END CLASS

vhf offers the right CAM solution for every dental laboratory and every practice. Because whether you are dry or wet processing – we have the perfect machine for you, with which you are optimally positioned for the future.



E3 E4 E5

N4+ K5 K5+ S5

Z4 R5

Material form

Disc



Blocks



Abutments



Thermo Forming Foil



Dry <sup>2</sup>



Wet <sup>2</sup>



Operation without  
ext. compressed air



Indication variety





Comfort



Typical field  
of application

Practice & practice lab



Lab & milling center



Detailed product  
information



<sup>1</sup> in preparation | <sup>2</sup> suitable materials for dry processing: plastics, wax, model plaster, zirconia, composites, CoCr • suitable materials for wet processing: glass ceramics, titanium • suitable materials for wet and dry processing: plastics, zirconia, composites, CoCr • the exact overview of the released materials per machines can be found in the respective factsheet.  
Be sure to review local and/or national regulations and/or regulations by other authorized organizations or entities (e.g. professional associations, health authorities).



## CREATING PERFECTION.

**vhf – synonymous with innovation and perfection since 1988.**

With over 35 years of experience in mechanical engineering, vhf is one of the leading manufacturers of dental milling machines. As a full-service CAM provider, vhf carefully develops and produces every single milling machine as well as the perfectly matched tools and software completely in-house. Everything from a single source. Made in Germany.

**Service. A matter close to our hearts.**

Despite their short maintenance intervals and particularly long service lives, servicing your machines is very important to us. We support you with our user-friendly dental **portal**, numerous online tutorials and personal support through our international service network.

### **vhf camufacture AG**

Lettenstraße 10  
72119 Ammerbuch  
Germany  
+49 7032 97097 000  
info@vhf.de | vhf.com

### **North America**

vhf Inc.  
80 Davids Drive, Suite 5  
Hauppauge, NY 11788, USA  
+1 631 524 5252  
info@vhf.com | vhf.com

### **Asia**

vhf Trading (Shanghai) Co., Ltd.  
Room 2902, Building T1, Tianshan SOHO,  
No. 421 Ziyun Road, Changning District,  
Shanghai, China  
asia@vhf.de | vhf.com

EN



EN  
DE  
FR  
IT  
ES  
CN

As of: 03/2025  
No. 269748

**vhf**  
CREATING PERFECTION