

Bestseller mills at the Greater New York Dental Meeting

vhf presents broad product portfolio at the largest U.S. dental show

Ammerbuch/Hauppauge, Nov. 23, 2018: vhf will present its current bestsellers during this year's Greater New York Dental Meeting from November 25–28, 2018. With the R5, Z4 and K5+, vhf is showcasing three different machines for wet and dry processing. The company is thereby further consolidating its position in the U.S. dental market, which has recently shown strong demand for vhf milling and grinding machines.

vhf will give an overview of its current four- and five-axis milling and grinding machines for dental labs and offices at booth 1317 of the German Pavilion in New York. The company offers outstanding solutions for both dry and wet processing. According to the vhf basic principle of open systems, users have free choice regarding components such as scanners and design software in addition to the materials to be processed.

With the futuristic Z4, vhf is showcasing a milling and grinding machine optimized for same-day restorations that fits perfectly into a digital workflow. Thanks to its ease-of-use, fast spindle and optimized machining strategies, restorations of the highest quality can be produced with the machine in minutes. This makes it possible to provide patients with durable, high-quality dental prostheses in a single day.

The R5 can be equipped with every type of material commonly used in dental technology in both disc and block form and up to a thickness of 40 mm. The highly automated five-axis milling and grinding machine for wet and dry processing comes with a variety of innovative features such as a patent-pending system for the quick tool-free loading of up to ten discs or 60 blocks or abutments as well as an automatic cleaning and drying system of the working chamber.

The latest milling machine from vhf is also represented in the form of the K5+, a new edition of the bestselling K5. The device is equipped with many new user-friendly features, such as tool-free blank loading (DirectDiscTechnology). Moreover, an integrated ionizer neutralizes the static charge of acrylic chips, and an integrated webcam simplifies support.

Images



The three vhf dental milling machines exhibited at GNYDM: K5+, Z4, and R5 (from left).

About vhf:

Founded in 1988, vhf camfacture AG is a leading manufacturer of CNC milling machines, milling tools and CAM software for the dental sector, industry and sign making. Headquartered in Ammerbuch, 30 kilometers southwest of Stuttgart, vhf employs more than 220 people and is constantly expanding. With its subsidiary in the state of New York, it offers services to U.S. customers.

Global:

vhf camfacture AG · Lettenstrasse 10 · D-72119 Ammerbuch · Germany
+49 7032 97097 000 · info@vhf.de · vhf.de

The Americas:

vhf Inc. · 80 Davids Drive, Suite 5 · Hauppauge, NY 11788 · USA
+1 631 524 5252 · info@vhf.com · vhf.com

About the Greater New York Dental Meeting:

The Greater New York Dental Meeting (GNYDM) is one of the largest trade shows in the industry worldwide and is taking place for the 94th time. Most recently, the show attracted more than 70,000 visitors, including almost 20,000 dentists and 9,000 international guests from 150 countries.

###

Download

The text of this press release and an overview of all vhf press releases and images of the past years can be found at: news.vhf.de

Press Contact

Matthias Brehmke, B. A.
PR Manager

+49 7032 97097 608
matthias.brehmke@vhf.de

Christine McClymont
Dipl.-Medienwirtin (FH)
Head of Marketing and Communications
+49 7032 97097 606
christine.mcclymont@vhf.de

For immediate release.